

Probe HER2/TOP2A/CEN17 TriColor FISH Probe(V.2)
Product No. GC247

HER2

Probe Size Approximately 510kb

Fluorophore FITC

Location 17q12

Also known as NEU; NGL; ERBB2; TKR1; CD340; HER-2; VSCN2; MLN 19; MLN-19; c-ERB2; c-ERB-2; HER-2/neu; p185(erbB2)

TOP2A

Probe Size Approximately 450kb

Fluorophore TAMRA(Orange)

Location 17q21

Also known as TOP2; TP2A; TOPIIA; TOP2alpha

CEN17(D17Z1)

Fluorophore DEAC

Location 17p11.1-17q11.1

