




Probe EGFR/CEN7q Dual Color FISH Probe
Probe No. GC005

EGFR
Probe Size Approximately 320kb
Fluorophore TexRed
Location 7p12
Also known as ERBB; HER1; mENA; ERBB1; PIG61; NISBD2

CEN7q
Probe Size Approximately 800kb
Fluorophore FITC
Location 7q11.21



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.