

Probe MLL/CEN11p Dual Color FISH Probe
Product No. GC016

MLL

Probe Size Approximately 380kb

Fluorophore TexRed

Location 11q23.3




Also known as HRX; MLL; MLL1; TRX1; ALL-1; CXXC7; HTRX1; MLL1A; WDSTS; MLL-AF9; MLL/GAS7; TET1-MLL

CEN11p

Probe Size Approximately 630kb

Fluorophore FITC

Location 11p11.12



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.