

Probe PTEN/CEN10p Dual Color FISH Probe  
Product No. GC019


**PTEN**

Probe Size Approximately 470kb  
Fluorophore TexRed  
Location 10q23  
Also known as BZS; DEC; CWS1; GLM2; MHAM; TEP1; MMAC1; PTEN1; 10q23del




**CEN10p**

Probe Size Approximately 640kb  
Fluorophore FITC  
Location 10p11.22

PTEN



CEN10p



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.