

Probe APC/CEN5q Dual Color FISH Probe
Product No. GC021


APC

Probe Size Approximately 420kb
Fluorophore TexRed
Location 5q22.2
Also known as GS; DP2; DP3; BTPS2; DP2.5; PPP1R46




CEN5q

Probe Size Approximately 700kb
Fluorophore FITC
Location 5q11.2

APC



CEN5q



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.