

Probe HER3/CEN12p Dual Color FISH Probe
Product No. GC023

HER3

Probe Size Approximately 320kb

Fluorophore TexRed

Location 12q13



Also known as ERBB3; LCCS2; ErbB-3; c-erbB3; erbB3-S; MDA-BF-1; c-erbB-3; p180-ErbB3; p45-sErbB3; p85-sErbB3

CEN12p

Probe Size Approximately 530kb

Fluorophore FITC

Location 12p11.22



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.