


Probe UBE3A/CEN15q Dual Color FISH Probe
Product No. GC025



UBE3A
Probe Size Approximately 370kb
Fluorophore TexRed
Location 15q11.2
Also known as AS; ANCR; E6-AP; HPVE6A; EPVE6AP

CEN15q
Probe Size Approximately 680kb
Fluorophore FITC
Location 15q13.3

UBE3A



CEN15q



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.