




Probe TRF1/CEN8p Dual Color FISH Probe
Product No. GC028

TRF1

Probe Size Approximately 290kb
Fluorophore TexRed
Location 8q21.11
Also known as TRF; PIN2; TRF1; TRBF1; t-TRF1; hTRF1-AS

CEN8p

Probe Size Approximately 520kb
Fluorophore FITC
Location 8p11.21



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.