

Probe VHL/CEN3q Dual Color FISH Probe
Product No. GC032


VHL

Probe Size Approximately 170kb
Fluorophore TexRed
Location 3p26-p25
Also known as RCA1; VHL1; pVHL; HRCA1




CEN3q

Probe Size Approximately 500kb
Fluorophore FITC
Location 3q12.1

VHL



CEN3q



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.