

Probe CALR/CEN19q Dual Color FISH Probe
Product No. GC033

CALR

Probe Size Approximately 260kb
Fluorophore TexRed
Location 19p13.1-13.2
Also known as RO; CRT; SSA; cC1qR; HEL-S-99n



CEN19q

Probe Size Approximately 430kb
Fluorophore FITC
Location 19q12

CALR



CEN19q



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.