

Probe CDK4/CEN12q Dual Color FISH Probe
Product No. GC038


CDK4

Probe Size Approximately 370kb
Fluorophore TexRed
Location 12q14.1
Also known as CMM3; PSK-J3

CEN12q

Probe Size Approximately 580kb
Fluorophore FITC
Location 12q12

CDK4





~370kb

CEN12q



~580kb



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.