


Probe YES1/CEN18q FISH Probe
 Product No. GC040




YES1
 Probe Size Approximately 150kb
 Fluorophore TexRed
 Location 18p11.32
 Also known as Yes; c-yes; HsT441; P61-YES

CEN18q
 Probe Size Approximately 660kb
 Fluorophore FITC
 Location 18q12.1

YES1



CEN18q



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.