

Probe TPX2/CEN20p Dual Color FISH Probe  
Product No. GC047

**TPX2**

Probe Size Approximately 150kb

Fluorophore TexRed

Location 20q11.2




Also known as DIL2; p100; DIL-2; HCTP4; FLS353; HCA519; REPP86; C20orf1; C20orf2; GD:C20orf1

**CEN20p**

Probe Size Approximately 400kb

Fluorophore FITC

Location 20p11.21



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.