

Probe DTNB/CEN2p Dual Color FISH Probe  
Product No. GC051

**DTNB**


Probe Size Approximately 310kb  
Fluorophore TexRed  
Location 2p24

**CEN2p**

Probe Size Approximately 670kb  
Fluorophore FITC  
Location 2p11.2



DTNB

25.2 25.4 25.6 25.8 26.0 26.2 26.4 (MB)



CEN2p

85.0 85.2 85.4 85.6 85.8 86.0 86.2 (MB)



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.