




Probe MYCN/CEN2p Dual Color FISH Probe
Product No. GC054

MYCN

Probe Size Approximately 270kb
Fluorophore TexRed
Location 2p24.1
Also known as NMYC; ODED; MODED; N-myc; bHLHe37

CEN2p

Probe Size Approximately 670kb
Fluorophore FITC
Location 2p11.2



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.