

Probe MDM4/CEN1q Dual Color FISH Probe
Product No. GC059


MDM4

Probe Size Approximately 350kb
Fluorophore TexRed
Location 1q32
Also known as HDMX; MDMX; MRP1




CEN1q

Probe Size Approximately 600kb
Fluorophore FITC
Location 1q21.3

MDM4



CEN1q



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.