

Probe ATP7B/CEN13q Dual Color FISH Probe
Product No. GC062


ATP7B

Probe Size Approximately 350kb
Fluorophore TexRed
Location 13q14.3
Also known as WD; PWD; WC1; WND




CEN13q

Probe Size Approximately 550kb
Fluorophore FITC
Location 13q12.11

ATP7B



CEN13q



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.