

Probe AKT3/CEN1q Dual Color FISH Probe
Product No. GC071

AKT3

Probe Size Approximately 600kb

Fluorophore TexRed

Location 1q43-q44

Also known as MPPH; PKBG; MPPH2; PRKBG; STK-2; PKB-GAMMA; RAC-gamma;
RAC-PK-gamma


CEN1q

Probe Size Approximately 600kb




Fluorophore FITC

Location 1q21.3

AKT3



CEN1q



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.