

Probe BRCA1/CEN17q Dual Color FISH Probe
Product No. GC073

BRCA1

Probe Size Approximately 150kb

Fluorophore TexRed

Location 17q21.31




Also known as IRIS; PSCP; BRCAI; BRCC1; FANCS; PNCA4; RNF53; BROVCA1; PPP1R53

CEN17q

Probe Size Approximately 540kb

Fluorophore FITC

Location 17q12



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.