

Probe CDH1/CEN16q Dual Color FISH Probe  
Product No. GC074


**CDH1**

Probe Size Approximately 390kb  
Fluorophore TexRed  
Location 16q22.1  
Also known as UVO; CDHE; ECAD; LCAM; Arc-1; CD324




**CEN16q**

Probe Size Approximately 700kb  
Fluorophore FITC  
Location 16q12.1

**CDH1**



**CEN16q**



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.