

Probe SOX2/CEN3q Dual Color FISH Probe
Product No. GC078


SOX2

Probe Size Approximately 380kb
Fluorophore TexRed
Location 3q26.3-q27
Also known as ANOP3; MCOPS3




CEN3q

Probe Size Approximately 500kb
Fluorophore FITC
Location 3q12.1

SOX2



CEN3q



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.