

Probe PDGFRA/CEN4p Dual Color FISH Probe
Product No. GC091


PDGFRA

Probe Size Approximately 240kb
Fluorophore TexRed
Location 4q12
Also known as CD140A; PDGFR2; PDGFR-2; RHEPDGFRA




CEN4p

Probe Size Approximately 560kb
Fluorophore FITC
Location 4p12

PDGFRA



CEN4p



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.