

Probe MER/CEN2q Dual Color FISH Probe
Product No. GC094


MER

Probe Size Approximately 290kb
Fluorophore TexRed
Location 2q13
Also known as MER; RP38; c-Eyk; c-mer; Tyro12




CEN2q

Probe Size Approximately 580kb
Fluorophore FITC
Location 2q11.2

MER



CEN2q



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.