

Probe BAP1/CEN3q Dual Color FISH Probe
Product No. GC097

BAP1

Probe Size Approximately 150kb
Fluorophore TexRed
Location 3p21.1
Also known as UCHL2; hucep-6; HUCEP-13

CEN3q


Probe Size Approximately 500kb
Fluorophore FITC
Location 3q12.1

BAP1





~150kb

CEN3q



~500kb



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.