

Probe MST1R/CEN3q Dual Color FISH Probe
Product No. GC098


MST1R

Probe Size Approximately 150kb
Fluorophore TexRed
Location 3p21.3
Also known as RON; PTK8; CD136; CDw136



CEN3q

Probe Size Approximately 500kb
Fluorophore FITC
Location 3q12.1

MST1R



CEN3q



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.