


Probe
Product No. HOXA9/CEN7p Dual Color FISH Probe
GC109




HOXA9
Probe Size Approximately 190kb
Fluorophore TexRed
Location 7p15.2
Also known asHOX1; ABD-B; HOX1G; HOX1.7

CEN7p
Probe Size Approximately 700kb
Fluorophore FITC
Location 7p12.1

HOXA9



CEN7p



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.