

Probe LAMP2/CENXp Dual Color FISH Probe
Product No. GC064

LAMP2

Probe Size Approximately 250kb
Fluorophore TexRed
Location Xq24
Also known as LAMPB; CD107b; LAMP-2; LGP110

CENXp

Probe Size Approximately 550kb
Fluorophore FITC
Location Xp11.22

LAMP2

119.0 119.2 119.4 119.6 119.8 120.0 120.2 (MB)





~250kb

CENXp

48.8 49.0 49.2 49.4 49.6 49.8 50.0 (MB)



~550kb



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.