

Probe TERC/CEN3q Dual Color FISH Probe
Product No. GC120

TERC

Probe Size Approximately 230kb
Fluorophore TexRed
Location 3q26
Also known as TR; hTR; TRC3; DKCA1; PFBMFT2; SCARNA19




CEN3q

Probe Size Approximately 500kb
Fluorophore FITC
Location 3q12.1

TERC



CEN3q



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.