

Probe TERT/CEN5q Dual Color FISH Probe
Product No. GC124


TERT

Probe Size Approximately 410kb
Fluorophore TexRed
Location 5p15.33
Also known as TP2; TRT; CMM9; EST2; TCS1; hTRT; DKCA2; DKCB4; hEST2; PFBMFT1




CEN5q

Probe Size Approximately 700kb
Fluorophore FITC
Location 5q11.2

TERT



CEN5q



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.