

Probe CDKN2C/CEN1q Dual Color FISH Probe
Product No. GC137


CDKN2C

Probe Size Approximately 150kb
Fluorophore TexRed
Location 1p32.3
Also known as p18; INK4C; p18-INK4C




CEN1q

Probe Size Approximately 600kb
Fluorophore FITC
Location 1q21.3

CDKN2C



CEN1q



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.