

Probe PRDM16/CEN1p Dual Color FISH Probe
Product No. GC142


PRDM16

Probe Size Approximately 450kb
Fluorophore TexRed
Location 1p36.23-33
Also known as MEL1; LVNC8; PFM13; CMD1LL




CEN1p

Probe Size Approximately 780kb
Fluorophore FITC
Location 1p13.3

PRDM16



CEN1p



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.