

Probe SRC/CEN20p Dual Color FISH Probe
Product No. GC149





SRC

Probe Size Approximately 150kb
Fluorophore TexRed
Location 20q12-q13
Also known as ASV; SRC1; c-SRC; p60-Src

CEN20p

Probe Size Approximately 400kb
Fluorophore FITC
Location 20p11.21

SRC



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.