

Probe WDFY2/CEN13q Dual Color FISH Probe
Product No. GC150

WDFY2

Probe Size Approximately 270kb
Fluorophore TexRed
Location 13q14.3
Also known as PROF; WDF2; ZFYVE22

CEN13q

Probe Size Approximately 550kb
Fluorophore FITC
Location 13q12.11

WDFY2

51.6 51.8 52.0 52.2 52.4 52.6 52.8 (MB)





~270kb

CEN13q

27.4 27.6 27.8 28.0 28.2 28.4 28.6 (MB)



~550kb



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.