

Probe GSTP1/CEN11p Dual Color FISH Probe
Product No. GC166





GSTP1

Probe Size Approximately 180kb
Fluorophore TexRed
Location 11q13.2
Also known as PI; DFN7; GST3; GSTP; FAEES3; HEL-S-22

CEN11p

Probe Size Approximately 630kb
Fluorophore FITC
Location 11p11.12

GSTP1



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.