

Probe SMARCB1/CEN22q Dual Color FISH Probe
Product No. GC181

SMARCB1

Probe Size Approximately 250kb

Fluorophore TexRed

Location 22q11.23

Also known as RDT; CSS3; INI1; SNF5; Snr1; BAF47; MRD15; RTPS1; Sfh1p; hSNFS; SNF5L1; SWNTS1; PPP1R144


CEN22q

Probe Size Approximately 620kb




Fluorophore FITC

Location 22q11.22

SMARCB1



CEN22q



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.