

Probe SP7/CEN12p Dual Color FISH Probe  
Product No. GC188


**SP7**

Probe Size Approximately 160kb  
Fluorophore TexRed  
Location 12q13.13  
Also known as OSX; OI11; OI12; osterix




**CEN12p**

Probe Size Approximately 530kb  
Fluorophore FITC  
Location 12p11.22

SP7



CEN12p



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.