

Probe MTAP/CEN9q Dual Color FISH Probe
Product No. GC206


MTAP

Probe Size Approximately 210kb
Fluorophore TexRed
Location 9p21.3
Also known as BDMF; MSAP; DMSFH; LGMBF; DMSMFH; c86fus; HEL-249




CEN9q

Probe Size Approximately 470kb
Fluorophore FITC
Location 9q21

MTAP



CEN9q



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.