





Probe HER2/CEN17 Dual Color FISH Probe
Product No. GC207

HER2


Probe Size Approximately 220kb
Fluorophore TexRed
Location 17q21.1
Also known as NEU; NGL; ERBB2; TKR1; CD340; HER-2; MLN 19; HER-2/neu

CEN17 (D17Z1)

Fluorophore FITC
Location 17p11.1-17q11.1



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.



Representative images of breast cancer stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.