

Probe TNFRSF8/CEN1p Dual Color FISH Probe
Product No. GC223


TNFRSF8

Probe Size Approximately 350kb
Fluorophore TexRed
Location 1p36.22
Also known as CD30; Ki-1; D1S166E




CEN1p

Probe Size Approximately 780kb
Fluorophore FITC
Location 1p13.3

TNFRSF8



CEN1p



representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.