

Probe BATF3/CEN1q Dual Color FISH Probe
Product No. GC224


BATF3

Probe Size Approximately 310kb
Fluorophore TexRed
Location 1q32.3
Also known as JDP1; SNFT; JUNDM1




CEN1q

Probe Size Approximately 600kb
Fluorophore FITC
Location 1q21.3

BATF3



CEN1q



representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.