

Probe CDK4/CEN12 Dual Color FISH Probe
Product No. GC232


CDK4

Probe Size Approximately 370kb
Fluorophore TexRed
Location 12q14.1
Also known as CMM3; PSK-J3




CEN12(D12Z3)

Probe Size Approximately 3Mb
Fluorophore FITC
Location 12p11.1-q11

CDK4



CEN12(D12Z3)



representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.