

Probe PLCG1/CEN20p Dual Color FISH Probe
Product No. GC234


PLCG1

Probe Size Approximately 300kb
Fluorophore TexRed
Location 20q12
Also known as PLC1; NCKAP3; PLC-II; PLC148; PLCgamma1



CEN20p

Probe Size Approximately 400kb
Fluorophore FITC
Location 20p11.21

PLCG1



CEN20p



representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.