

Probe VEGFR3/CEN5q Dual Color FISH Probe
Product No. GC237

VEGFR3

Probe Size Approximately 160kb
Fluorophore TexRed
Location 5q35.3
Also known as PCL; CHTD7; FLT-4; FLT41; LMPH1A; LMPHM1; VEGFR-3

CEN5q

Probe Size Approximately 700kb
Fluorophore FITC
Location 5q11.2

VEGFR3

179.4 179.6 179.8 180.0 180.2 180.4 180.6 (MB)





~160kb

CEN5q

51.6 51.8 52.0 52.2 52.4 52.6 52.8 (MB)



~700kb



representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.