

Probe HER2/CEN17 Dual Color FISH Probe (V.2)
Product No. GC241

HER2

Probe Size Approximately 510kb

Fluorophore Orange



Location 17q12

Also known as NEU; NGL; HER2; TKR1; CD340; HER-2; VSCN2; MLN 19; HER-2/neu

CEN17 (D17Z1)

Fluorophore FITC

Location 17p11.1-17q11.1



representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.