

Probe GAPDH/CEN12q FISH Probe
Product No. GC253


GAPDH

Probe Size Approximately 530kb
Fluorophore TexRed
Location 12p13.31
Also known as G3PD; GAPD; HEL-S-162eP




CEN12q

Probe Size Approximately 580kb
Fluorophore FITC
Location 12q12

GAPDH



CEN12q



representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.