

Probe WNT10B/CEN12p Dual Color FISH Probe
Product No. GC258


WNT10B

Probe Size Approximately 360kb
Fluorophore TexRed
Location 12q13.12
Also known as SHFM6; STHAG8; WNT-12




CEN12p

Probe Size Approximately 530kb
Fluorophore FITC
Location 12p11.22

WNT10B



CEN12p



representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus