

Probe TPSAB1/CEN16q Dual Color FISH Probe
Product No. GC263


TPSAB1

Probe Size Approximately 150kb
Fluorophore TexRed
Location 16p13.3
Also known as TPS1; TPS2; TPSB1; TPSB2; Tryptase-2




CEN16q

Probe Size Approximately 700kb
Fluorophore FITC
Location 16q12.1

TPSAB1



CEN16q



representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus