

Probe LAMC2/CEN1q Dual Color FISH Probe  
Product No. GC259


#### LAMC2

Probe Size Approximately 300kb  
Fluorophore TexRed  
Location 1q25.3  
Also known as B2T; CSF; EBR2; BM600; EBR2A; JEB3A; JEB3B; LAMB2T; LAMNB2




#### CEN1q

Probe Size Approximately 600kb  
Fluorophore FITC  
Location 1q21.3

#### LAMC2



#### CEN1q



representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus