

Probe RB1/CEN13q Dual Color FISH Probe
Product No. GC132

RB1

Probe Size Approximately 260kb
Fluorophore TexRed
Location 13q14.2
Also known as RB; pRb; OSRC; pp110; p105-Rb; PPP1R130

CEN13q

Probe Size Approximately 550kb
Fluorophore FITC
Location 13q12.11

RB1


48.4 48.6 48.8 49.0 49.2 49.4 49.6 (MB)





~260kb

CEN13q

27.4 27.6 27.8 28.0 28.2 28.4 28.6 (MB)



~550kb



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.