

Probe PPARG/CEN3q Dual Color FISH Probe
Product No. GC189


PPARG

Probe Size Approximately 190kb
Fluorophore TexRed
Location 3p25
Also known as GLM1; CIMT1; NR1C3; PPARG1; PPARG2; PPARGgamma




CEN3q

Probe Size Approximately 500kb
Fluorophore FITC
Location 3q12.1

PPARG



CEN3q



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.